

NPGS MARKER AND SPECULUMS

MERMOUD

- 4.6999** Non Penetrating
Glaucoma
Surgery Marker
4 x 5 mm
and 4 x 4 mm



BARRAQUER

- Eye speculum
- 4.7030** Baby, 3 mm
- 4.7000** Baby, 5 mm
- 4.7001** Child, 10 mm
- 4.7028** Adult, 12 mm
- 4.7002** Adult, 15 mm
- 4.7010** Adult, 15 mm
heavy wire



CORCELLE

- Eye speculum
- 4.7003** Baby, 8 mm
- 4.7004** Child, 10 mm
- 4.7005** Adult, 12 mm



CASTROVIEJO

- 4.7007** Eye speculum,
child
- 4.7008** Adult



BANGERTER


- 4.7014** Eye speculum,
adjustable spread



SHAARAWY

- 4.7013** Eye speculum W/4 nobs
permits for attachment of
corneal traction sutures





SPECULUMS - RETRACTORS

KNORZ

4.7019 Lakik eye speculum



GRAUWELS

4.7023 Horse speculum
Length: 19 cm
Spread: 11 cm
Blades: 3.4 cm




E = 1/5

KRATZ-BARRAQUER

Speculum

4.7025 Heavy wire, large
4.7026 Standard wire, large
4.7027 Standard wire, small




BARRAQUER

4.7029 Speculum
with solid blades,
large



LIEBERMANN

4.7035 Speculum KRATZ-Wire
Adjustable with 15.0 mm
open blades
Angled at 45 degrees



4.7036 With Aspiration

STEVENSON

4.7045 Lacrimal sac retractor
2 x 3, blunt prongs self
retaining adjustable spread




HARTMANN - HALLE

Nasal speculum

4.7047 Infant
4.7048 Child
4.7049 Adult



E = 1/2



LID PLATE - DILATOR - RETRACTORS

JAEGER

4.7100 Lid plate, stainless steel



SCHOCKET

4.7150 Scleral depressor, double end



DESMARRES

Lid retractor

- 4.7200** 10 mm
- 4.7201** 12 mm
- 4.7202** 14 mm
- 4.7203** 16 mm



BEEHLER

Pupil dilator 2.8 mm for 2.5 mm incision

- 4.7245** Straight
- 4.7246** Curved
- 4.7247** Replacement kit



KNAPP

- 4.7250** Lacrimal retractor, 3 prongs, 6 mm
- 4.7251** 4 prongs, 8 mm



ROTH

4.7260 Muscle retractor, 6 mm



BONN

4.7270 Orbital retractor 10 mm

