


STERILIZING AND STORAGE TRAYS UP TO 190°C

- Provides maximum protection as well as ease of fixation for a variety of instruments
- Complete instruments set are visible
- Ease of transportation to surgical sites
- Withstands autoclave and dry heat up to 190°C as well as EO and cold solutions

STERILISABLE PLASTIC TRAYS

- 4.9202** For 2 "HUCO" knives
with racks,
O.D.: 58 x 165 x
25 mm




- 4.9203** For 2-3 instruments
with mat,
O.D.: 84 x 165 x
25 mm



- 4.9210** For approximately
10 instruments
with mat,
O.D.: 178 x 285 x
30 mm







STERILIZING AND STORAGE TRAYS


STERILISABLE PLASTIC TRAYS

- 4.9220** For approximately 20 instruments, O.D.: 244 x 389 x 30 mm
- 4.9230** For approximately 30 instruments, O.D.: 335 x 536 x 30 mm




-
- 4.9301** For 1 large item, O.D. 149 x 249 x 40 mm
 - 4.9310** For approximately 10 instruments, with 1 silicone mat, O.D. 149 x 249 x 23 mm
 - 4.9320** For approximately 20 instruments, with 1 insert and 2 silicone mats, O.D. 149 x 249 x 40 mm
 - 4.9325** For approximately 10-15 instruments, with 8 bars to be cut by the user according to what instruments will be placed in the tray, O.D. 149 x 249 x 40 mm
 - 4.9332** Silicone mat only, 126 x 215 x 16 mm






STERILIZING AND STORAGE TRAYS

STERILISABLE PLASTIC TRAYS



- 4.9201** For 1-2 instruments with 1 silicone mat, O.D. 58 x 165 x 25 mm
- 4.9203** For 2-3 instruments with 1 silicone mat, O.D. 84 x 165 x 25 mm
- 4.9310** For approximately 10 instruments with 1 silicone mat, O.D. 149 x 250 x 25 mm
- 4.9331** Silicone mat only tray 4.9310
- 4.9312** For approximately 10 instruments with 1 insert and 1 silicone mat, O.D. 149 x 250 x 40 mm
- 4.9332** Silicone mat for lower level tray 4.9312
- 4.9333** Silicone mat for insert tray 4.9312
- 4.9341** For approximately 20 instruments with 1 silicone mat, O.D. 267 x 394 x 34 mm
- 4.9343** For approximately 25 instruments with 1 insert and 1 silicone mat, O.D. 267 x 394 x 50 mm
- 4.9330** Silicone mat for lower level tray 4.9343



STERILIZING AND STORAGE TRAYS

STERILISABLE PLASTIC TRAYS



- 4.9224 For 10 instruments with 1 silicone mat, O.D. 177 x 282 x 50 mm
- 4.9244 Silicone mat for lower level

- 4.9225 For 20 instruments with 1 silicone mat, O.D. 275 x 310 x 50 mm
- 4.9245 Silicone mat for lower level

- 4.9226 For 20 instruments with 1 silicone mat, O.D. 246 x 390 x 45 mm
- 4.9246 Silicone mat for lower level

- 4.9227 For 30 instruments with 1 silicone mat, O.D. 335 x 537 x 53 mm
- 4.9232 Silicone mat for lower level